



FreeStyle

Libre 2

PLUS



Useful tips while using FreeStyle Libre 2 Plus

- 1 Store the Sensor Pack in a temperature range from 4 °C to 25 °C in a cool dry area. You can also put the sensor in the refrigerator as long as the temperature is in the recommended range
- 2 Only apply to the back of your upper arm. Avoid areas with scars, moles, stretch marks, or lumps
- 3 You can scan your FreeStyle Libre 2 Plus sensor through clothing with a thickness of 4 cm
- 4 Wash the area with a non-moisturising, fragrance-free soap. Moisturising soaps can leave a film that prevents the adhesive from sticking

- 5 The sensor is water-resistant up to 3 feet (1 meter) for 30 minutes. Pat it dry gently with a towel after swimming or showering. Take extra care when towelling to avoid catching or pulling your sensor
- 6 Give your sensor room to breathe by wearing loose-fitting clothing.
- 7 Be careful not to catch your sensor on a doorway, car door, seat belt or furniture edges
- 8 Never use the exact same spot twice in a row. Move at least 1 inch from the previous site to prevent lipohypertrophy
- 9 Check your airline's policy on NFC and Bluetooth devices



Scan to visit our website
www.FreeStyleLibre.co.in

Got a question on your FreeStyle Libre 2 Plus system? We're here to help.



1800 2100 998



adc@cs.in@abbott.com

**Libre Customer Care Hours: 8:30am – 5:30pm,
Monday to Friday, excluding Saturday, Sunday and national public holidays.**

For the use of registered medical practitioners or hospitals or pharmacies. Images are for illustration purpose only. No actual patient data. Information contained herein is for distribution outside of the USA only.

*The FreeStyle LibreLink app is only compatible with certain mobile devices and operating systems. Please check the website for more information about device compatibility before using the app. Use of FreeStyle LibreLink may require registration with LibreView.

ADC Corporate Address: Abbott Healthcare Pvt. Ltd., Floor 16, Godrej BKC, Plot No. C- 68, BKC, Near MCA Club, Bandra (E), Mumbai-400 051. Registered Address: Abbott Healthcare Pvt. Ltd., 3, Corporate Park, Sion-Trombay Road, Chembur, Mumbai-400071.

© 2026 Abbott. All rights reserved. The sensor housing, FreeStyle, Libre, and related brand marks are marks of Abbott. ADC-2693099 v1.0 08/26

